

## Product datasheet for **RC214469L4V**

### ARTS1 (ERAP1) (NM\_016442) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	ARTS1 (ERAP1) (NM_016442) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ARTS1
Synonyms:	A-LAP; ALAP; APPILS; ARTS-1; ARTS1; ERAAP; ERAAP1; PILS-AP; PILSAP
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_016442
ORF Size:	2844 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC214469).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_016442.3</a>
RefSeq Size:	5584 bp
RefSeq ORF:	2847 bp
Locus ID:	51752
UniProt ID:	<a href="#">Q9NZ08</a>
Cytogenetics:	5q15
Domains:	Peptidase_M1
Protein Families:	Druggable Genome, Protease, Secreted Protein



[View online »](#)

**MW:** 107.7 kDa

**Gene Summary:** The protein encoded by this gene is an aminopeptidase involved in trimming HLA class I-binding precursors so that they can be presented on MHC class I molecules. The encoded protein acts as a monomer or as a heterodimer with ERAP2. This protein may also be involved in blood pressure regulation by inactivation of angiotensin II. Three transcript variants encoding two different isoforms have been found for this gene.[provided by RefSeq, Oct 2010]