

## Product datasheet for **RC214473L4V**

### CNBP (NM\_003418) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	CNBP (NM_003418) Human Tagged ORF Clone Lentiviral Particle
Symbol:	CNBP
Synonyms:	CNBP1; DM2; PROMM; RNF163; ZCCHC22; ZNF9
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_003418
ORF Size:	531 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC214473).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_003418.1</a>
RefSeq Size:	1500 bp
RefSeq ORF:	534 bp
Locus ID:	7555
UniProt ID:	<a href="#">P62633</a>
Cytogenetics:	3q21.3
Domains:	zf-CCHC
Protein Families:	Druggable Genome, Transcription Factors



[View online »](#)

**MW:** 19.3 kDa

**Gene Summary:** This gene encodes a nucleic-acid binding protein with seven zinc-finger domains. The protein has a preference for binding single stranded DNA and RNA. The protein functions in cap-independent translation of ornithine decarboxylase mRNA, and may also function in sterol-mediated transcriptional regulation. A CCTG expansion from <30 repeats to 75-11000 repeats in the first intron of this gene results in myotonic dystrophy type 2. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2016]