

Product datasheet for **RC214501A1V**

Human Epigen (EPGN) (NM_001013442) AAV Particle

Product data:

Product Type: AAV Particles

Product Name: Human Epigen (EPGN) (NM_001013442) AAV Particle

Tag: Myc-DDK

Symbol: Epigen

Synonyms: ALGV3072; EPG; epigen; PRO9904

Mammalian Cell Selection: None

Vector: pAAV-AC-Myc-DDK (PS100089)

ORF Nucleotide Sequence: >RC214501 representing NM_001013442
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGCTTTGGGAGTTCCAATATCAGTCTATCTTTTATTCAACGCAATGACAGCACTGACCGAAGAGGCAG
CCGTGACTGTAACACCTCCAATCACAGCCAGCAAGCTGACAACATAGAAGACCCATAGCCTTGAAGTT
CTCACACCTTTGCCTGGAAGATCATAACAGTTACTGCATCAACGGTGCTTGTGCATTCCACCATGAGCTA
GAGAAAGCCATCTGCAGGTGTTTTACTGGTTTACTGGAGAAAGGTGTGAGCACTTGACTTTAATTTCAT
ATGCTGTGGATTCTTATGAAAAATACATTGCAATTGGGATTGGTGTGGATTACTATTAAGTGGTTTTCT
TGTTATTTTTTACTGCTATATAAGAAAGAGGTATGAAAAAGACAAAATA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC214501 representing NM_001013442
Red=Cloning site Green=Tags(s)

MALGVPIISVYLLFNAMTALTEEAATVTPPITAQQADNIEGPALKFSLCLLEDHNSYINGACAFHHEL
EKAICRCFTGYTGERCEHLTLTSYAVDSYEKYIAIGIGVGLLLSGFLVIFYCYIRKRYEKDKI

TRTRPLE**QKL**ISEEDLAANDILDYKDDDDKV

Species: Human

Serotype: AAV-2

ACCN: NM_001013442



[View online »](#)

ORF Size:	399 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001013442.1, NP_001013460.1</u>
RefSeq Size:	847 bp
RefSeq ORF:	401 bp
Locus ID:	255324
Cytogenetics:	4q13.3
MW:	14.6 kDa