

## Product datasheet for **RC214517L4V**

### **NRN1L (NM\_198443) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	NRN1L (NM_198443) Human Tagged ORF Clone Lentiviral Particle
Symbol:	NRN1L
Synonyms:	cpg15-2; MRCC2446; UNQ2446
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_198443
ORF Size:	495 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC214517).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_198443.1</a> , <a href="#">NP_940845.1</a>
RefSeq Size:	633 bp
RefSeq ORF:	498 bp
Locus ID:	123904
UniProt ID:	<a href="#">Q496H8</a>
Cytogenetics:	16q22.1
MW:	17.8 kDa



[View online »](#)

**Gene Summary:**

The protein encoded by this gene is extracellular and enhances both neurite growth and neuronal survival. The encoded protein is found both as a GPI anchored membrane-bound form and as a secreted form. This activity-related ligand functions as a homodimer or heterodimer. [provided by RefSeq, Feb 2017]