

Product datasheet for **RC214539L4V**

PIGG (NM_017733) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	PIGG (NM_017733) Human Tagged ORF Clone Lentiviral Particle
Symbol:	PIGG
Synonyms:	GPI7; LAS21; MRT53; PRO4405; RLGS1930
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_017733
ORF Size:	2925 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC214539).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_017733.3
RefSeq Size:	3242 bp
RefSeq ORF:	2928 bp
Locus ID:	54872
UniProt ID:	Q5H8A4
Cytogenetics:	4p16.3
Protein Families:	Transmembrane
Protein Pathways:	Glycosylphosphatidylinositol(GPI)-anchor biosynthesis



[View online »](#)

MW: 107.4 kDa

Gene Summary: This gene encodes an enzyme involved in glycosylphosphatidylinositol-anchor biosynthesis. The encoded protein, which is localized to the endoplasmic reticulum, is involved in transferring ethanolamine phosphate to mannose 2 of glycosylphosphatidylinositol species H7 to form species H8. Allelic variants of this gene have been associated with intellectual disability, hypotonia, and early-onset seizures. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Sep 2016]