

Product datasheet for RC214541

CCL14 (NM_032962) Human Tagged ORF Clone

Product data:

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn


Product Type:	Expression Plasmids
Product Name:	CCL14 (NM_032962) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	CCL14
Synonyms:	CC-1; CC-3; CKB1; HCC-1; HCC-1(1-74); HCC-1/HCC-3; HCC-3; MCIF; NCC-2; NCC2; SCYA14; SCYL2; SY14
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
ORF Nucleotide Sequence:	>RC214541 representing NM_032962 Red=Cloning site Blue=ORF Green=Tags(s)
	TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC GCC <mark>GCGATCGC</mark> C
	ATGAAGATCTCCGTGGCTGCCATTCCCTTCTTCCTCCTCATCACCATCGCCCTAGGGACCAAGACTGAAT CCTCCTCACAAACTGGGGGGGAAACCGAAGGTTGTTAAAATACAGCTAAAGTTGGTGGGGGGGACCTTACCA CCCCTCAGAGTGCTGCTTCACCTACACTACCTACAAGATCCCGCGTCAGCGGATTATGGATTACTATGAG ACCAACAGCCAGTGCTCCAAGCCCGGAATTGTCTTCATCACCAAAAGGGGCCATTCCGTCTGTACCAACC CCAGTGACAAGTGGGTCCAGGACTATATCAAGGACATGAAGGAGAAC
	ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT ACAAGGATGACGACGATAAG GTTTAA
Protein Sequence:	>RC214541 representing NM_032962 Red=Cloning site Green=Tags(s)
	MKISVAAIPFFLLITIALGTKTESSSQTGGKPKVVKIQLKLVGGPYHPSECCFTYTTYKIPRQRIMDYYE TNSQCSKPGIVFITKRGHSVCTNPSDKWVQDYIKDMKEN
	TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Restriction Sites:	Sgfl-Mlul



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2022 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US




Cloning Scheme:



* The last codon before the Stop codon of the ORF

Plasmid Map:



ACCN:	
ORF Size:	
OTI Disclaimer:	

ACCNI

NM_032962

327 bp

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <u>More info</u>

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2022 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

CCL14	(NM_032962) Human Tagged ORF Clone – RC214541
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	 Centrifuge at 5,000xg for 5min. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. Close the tube and incubate for 10 minutes at room temperature. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM 032962.4, NP 116738.1</u>
RefSeq Size:	579 bp
RefSeq ORF:	330 bp
Locus ID:	6358
UniProt ID:	<u>Q16627</u>
Cytogenetics:	17q12
Protein Families:	Druggable Genome, Secreted Protein
Protein Pathways:	Chemokine signaling pathway, Cytokine-cytokine receptor interaction
MW:	12.3 kDa
Gene Summary:	This gene, chemokine (C-C motif) ligand 14, is one of several CC cytokine genes clustered on 17q11.2. The CC cytokines are secreted proteins characterized by two adjacent cysteines. The cytokine encoded by this gene induces changes in intracellular calcium concentration and enzyme release in monocytes. Multiple transcript variants encoding different isoforms have been found for this gene. Read-through transcripts are also expressed that include exons from the upstream cytokine gene, chemokine (C-C motif) ligand 15, and are represented as GeneID: 348249. [provided by RefSeq, Dec 2009]

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2022 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US