

Product datasheet for **RC214542L3V**

RGS11 (NM_003834) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	RGS11 (NM_003834) Human Tagged ORF Clone Lentiviral Particle
Symbol:	RGS11
Synonyms:	RS11
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_003834
ORF Size:	1338 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC214542).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_003834.1
RefSeq Size:	2384 bp
RefSeq ORF:	1341 bp
Locus ID:	8786
UniProt ID:	O94810
Cytogenetics:	16p13.3
Protein Families:	Druggable Genome
MW:	50.5 kDa


[View online »](#)

Gene Summary:

The protein encoded by this gene belongs to the RGS (regulator of G protein signaling) family. Members of the RGS family act as GTPase-activating proteins on the alpha subunits of heterotrimeric, signal-transducing G proteins. This protein inhibits signal transduction by increasing the GTPase activity of G protein alpha subunits, thereby driving them into their inactive GDP-bound form. Alternative splicing occurs at this locus and four transcript variants encoding distinct isoforms have been identified. [provided by RefSeq, Nov 2013]