

Product datasheet for RC214551

Product datasneet for RC214551

DEFB126 (NM_030931) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids

Product Name: DEFB126 (NM_030931) Human Tagged ORF Clone

Tag: Myc-DDK Symbol: DEFB126

Synonyms: bA530N10.1; C20orf8; DEFB-26; DEFB26; hBD-26; HBD26

Mammalian Cell N

Selection:

Neomycin

Vector:pCMV6-Entry (PS100001)E. coli Selection:Kanamycin (25 ug/mL)

ORF Nucleotide >RC214551 representing NM_030931

Sequence: Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC

GCCGCGATCGCC

ATGAAGTCCCTACTGTTCACCCTTGCAGTTTTTATGCTCCTGGCCCAATTGGTCTCAGGTAATTGGTATG
TGAAAAAGTGTCTAAACGACGTTGGAATTTGCAAGAAGAAGTGCAAACCTGAAGAGATGCATGTAAAGAA
TGGTTGGGCAATGTGCGGCAAACAAAGGGACTGCTGTTTCCAGCTGACAGACGTGCTAATTATCCTGTT
TTCTGTGTCCAGACAAAGACTACAAGAATTTCAACAGTAACAGCAACAACAACAACAACTTTGATGA

TGACTACTGCTTCGATGTCTTCGATGGCTCCTACCCCCGTTTCTCCCACTGGT

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT

ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC214551 representing NM_030931

Red=Cloning site Green=Tags(s)

MKSLLFTLAVFMLLAQLVSGNWYVKKCLNDVGICKKKCKPEEMHVKNGWAMCGKQRDCCVPADRRANYPV

 ${\tt FCVQTKTTRISTVTATTATTTLMMTTASMSSMAPTPVSPTG}$

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mk6112 b09.zip

Restriction Sites: Sgfl-Mlul



OriGene Technologies, Inc. 9620 Medical Center Drive, Ste 200


CN: techsupport@origene.cn

Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com



Cloning Scheme:





^{*} The last codon before the Stop codon of the ORF

ACCN: NM_030931

ORF Size: 333 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of

reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube

containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method: 1. Centrifuge at 5,000xg for 5min.

2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.

3. Close the tube and incubate for 10 minutes at room temperature.

4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid

at the bottom.

5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of

shipping when stored at -20°C.

RefSeq: NM 030931.4

RefSeq Size: 413 bp RefSeq ORF: 336 bp Locus ID: 81623



Cytogenetics:

UniProt ID: Q9BYW3

Protein Families: Protease, Secreted Protein, Transmembrane

20p13


MW: 10 kDa

Gene Summary: Defensins are cysteine-rich cationic polypeptides that are important in the immunologic


response to invading microorganisms. The antimicrobial protein encoded by this gene is secreted and is a member of the beta defensin protein family. Beta defensin genes are found in several clusters throughout the genome, with this gene mapping to a cluster at 20p13. The encoded protein is highly similar to an epididymal-specific secretory protein (ESP13.2) from cynomolgus monkey. Mutation of this gene is associated with impaired sperm function.

[provided by RefSeq, Nov 2014]

Product images:



Circular map for RC214551



Western blot validation of overexpression lysate (Cat# [LY403096]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from untransfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC214551 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).