

Product datasheet for **RC214566L3V**

SETD3 (NM_032233) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	SETD3 (NM_032233) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SETD3
Synonyms:	C14orf154; hSETD3
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_032233
ORF Size:	1782 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC214566).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_032233.2
RefSeq Size:	2905 bp
RefSeq ORF:	1785 bp
Locus ID:	84193
UniProt ID:	Q86TU7
Cytogenetics:	14q32.2
Domains:	SET
MW:	67.1 kDa



[View online »](#)

Gene Summary:

Protein-histidine N-methyltransferase that specifically mediates methylation of actin at 'His-73' (PubMed:30526847, PubMed:30626964, PubMed:30785395). Histidine methylation of actin is required for smooth muscle contraction of the laboring uterus during delivery (PubMed:30626964). Does not have protein-lysine N-methyltransferase activity and probably only catalyzes histidine methylation of actin (PubMed:30626964, PubMed:30785395). [UniProtKB/Swiss-Prot Function]