

## Product datasheet for **RC214579L1V**

### **MCH2 (MCHR2) (NM\_032503) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	MCH2 (MCHR2) (NM_032503) Human Tagged ORF Clone Lentiviral Particle
Symbol:	MCHR2
Synonyms:	GPR145; GPRv17; MCH-2R; MCH-R2; MCH2; MCH2R; MCHR-2; SLT
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_032503
ORF Size:	1020 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC214579).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_032503.2</a> , <a href="#">NP_115892.2</a>
RefSeq Size:	2144 bp
RefSeq ORF:	1023 bp
Locus ID:	84539
UniProt ID:	<a href="#">Q969V1</a>
Cytogenetics:	6q16.2
Domains:	7tm_1
Protein Families:	Druggable Genome, GPCR, Transmembrane



[View online »](#)

**Protein Pathways:** Neuroactive ligand-receptor interaction

**MW:** 38.8 kDa

**Gene Summary:** Receptor for melanin-concentrating hormone, coupled to G proteins that activate phosphoinositide hydrolysis.[UniProtKB/Swiss-Prot Function]