

Product datasheet for **RC214591L3V**

MGAT5B (NM_144677) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	MGAT5B (NM_144677) Human Tagged ORF Clone Lentiviral Particle
Symbol:	MGAT5B
Synonyms:	GnT-IX; GnT-VB
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_144677
ORF Size:	2370 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC214591).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_144677.2
RefSeq Size:	4246 bp
RefSeq ORF:	2373 bp
Locus ID:	146664
UniProt ID:	Q3V5L5
Cytogenetics:	17q25.2
Protein Families:	Transmembrane
Protein Pathways:	Metabolic pathways, N-Glycan biosynthesis



[View online »](#)

MW: 89.1 kDa

Gene Summary: The MGAT5B gene encodes a beta-1,6-N-acetylglucosaminyltransferase (EC 2.4.1.155) that functions in the synthesis of complex cell surface N-glycans (Kaneko et al., 2003 [PubMed 14623122]).[supplied by OMIM, Nov 2008]