

## Product datasheet for **RC214623L1V**

### **BAZ1A (NM\_013448) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	BAZ1A (NM_013448) Human Tagged ORF Clone Lentiviral Particle
Symbol:	BAZ1A
Synonyms:	ACF1; hACF1; WALp1; WCRF180
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_013448
ORF Size:	4668 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC214623).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_013448.2</a>
RefSeq Size:	6031 bp
RefSeq ORF:	4671 bp
Locus ID:	11177
UniProt ID:	<a href="#">Q9NRL2</a>
Cytogenetics:	14q13.1-q13.2
Protein Families:	Druggable Genome
MW:	178.5 kDa



[View online »](#)

**Gene Summary:**

The BAZ1A gene encodes the accessory subunit of the ATP-dependent chromatin assembly factor (ACF), a member of the ISWI ('imitation switch') family of chromatin remodeling complexes (summarized by Racki et al., 2009 [PubMed 20033039]).[supplied by OMIM, Apr 2010]