

Product datasheet for **RC214649L3V**

CCDC85A (NM_001080433) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	CCDC85A (NM_001080433) Human Tagged ORF Clone Lentiviral Particle
Symbol:	CCDC85A
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001080433
ORF Size:	1659 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC214649).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001080433.1
RefSeq Size:	3983 bp
RefSeq ORF:	1662 bp
Locus ID:	114800
UniProt ID:	Q96PX6
Cytogenetics:	2p16.1
MW:	59.8 kDa
Gene Summary:	May play a role in cell-cell adhesion and epithelium development through its interaction with proteins of the beta-catenin family.[UniProtKB/Swiss-Prot Function]



[View online »](#)