

## Product datasheet for **RC214676L3V**

### Stella (DPPA3) (NM\_199286) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Stella (DPPA3) (NM_199286) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Stella
Synonyms:	Pgc7; STELLA
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_199286
ORF Size:	477 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC214676).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_199286.2</a>
RefSeq Size:	1040 bp
RefSeq ORF:	480 bp
Locus ID:	359787
UniProt ID:	<a href="#">Q6W0C5</a>
Cytogenetics:	12p13.31
MW:	17.9 kDa


[View online »](#)

**Gene Summary:**

This gene encodes a protein that in mice may function as a maternal factor during the preimplantation stage of development. In mice, this gene may play a role in transcriptional repression, cell division, and maintenance of cell pluripotentiality. In humans, related intronless loci are located on chromosomes 14 and X. [provided by RefSeq, Jul 2008]