

Product datasheet for **RC214684A1V**

Human RING finger protein unkempt like (UNKL) (NM_001037125) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human RING finger protein unkempt like (UNKL) (NM_001037125) AAV Particle
Tag:	Myc-DDK
Symbol:	RING finger protein unkempt like
Synonyms:	C16orf28; ZC3H5L; ZC3HDC5L
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC214684 representing NM_001037125 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGCCGTCGGTTTCGAAAGCGGCGGCAGCGCGCTGAGCGGGTCCCCCGCAGACTGAGAAGCCGACCC
ACTACAGGTACCTGAAGGAGTTCAGGACGGAGCAGTGCCCCCTGTTTTCACAGCACAAGTGCGCGCAGCA
CCGGCCGTTACCTGCTTCCACTGGCACTTCTCAACCAGCGGCGCCGAGGCCCTCCGAGGCGCGAC
GGCACCTTCAACTACAGCCCCGACGTGTACTGCTCCAAGTACAACGAAGCCACCGCGGTGTGCCCGACG
GCGACGAGTGTCCCTACCTGCACCGGACGACGGGGACACAGAACGAAGTACCACCTGCGTTACTACAA
AACAGGAACCTGCATCCACGAGACAGACGACGTGGCCACTGCGTGAAGAATGGGCTGCACTGTGCCTTC
GCGCACGGCCCCCTGGACCTGCGGCCGCGCCGTGTGTGACGTGAGGGAGCTGACGGCCAGGAAGCCTTGC
AGAACGGCCAGCTGGGCGGCGGGGAAGGGGTCCCGGATCTGCAGCCTGGGGTCTTGCCAGCCAGGCCAT
GATTGAGAAGATCCTGAGCGAGGACCCCGGTGGCAAGATGCCAACTTCGTGCTGGGACGTACAAGACG
GAGCAGTGCCCGAAGCCGCCACGCTGTGCCGCGAGGGCTATGCGTGCCCACTACCACAATAGCCGGG
ACAGGCGGCGCAACCCCGGCGGTTCCAGTACAGCTGGCAGCTGGGACGCGGGTCTTAGGCTGAGTCC
CAGGGCCAACAACCCAAGGTCGCCCTGCCAGGGTGACACAGGACCTTCTCCACCGCC

ACGCGTACGCGGCGGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC214684 representing NM_001037125 Red=Cloning site Green=Tags(s) MPSVSKAAAAALSGSPQTEKPTHYRYLKEFRTEQCPLFSQHKCAQHRPFTCFHWHFLNQRRRPLRRRD GTFNYS PDVYCSKYNEATGVCPDGDECPYLHRTTGDTERKYHLRYKGTGCIHETDARGHCVKNGLHCAF AHGPLDLRPPVCDVRELQAQEALQNGQLGGGEGVPDLQPGVLASQAMIEKILSEDPWQDANFVLGSYKT EQCPKPPRLCRQGYACPHYHNSRDRRRNPRRFQYSWQLGRRVLRLSPRANNPRVALPRVHTGPSSTA TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_001037125
ORF Size:	831 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001037125.1</u>
RefSeq Size:	1401 bp
RefSeq ORF:	834 bp
Locus ID:	64718
UniProt ID:	<u>Q9H9P5</u>
Cytogenetics:	16p13.3
MW:	31.6 kDa