

## Product datasheet for **RC214684L4V**

### **RING finger protein unkempt like (UNKL) (NM\_001037125) Human Tagged ORF Clone Lentiviral Particle**

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	RING finger protein unkempt like (UNKL) (NM_001037125) Human Tagged ORF Clone Lentiviral Particle
Symbol:	RING finger protein unkempt like
Synonyms:	C16orf28; ZC3H5L; ZC3HDC5L
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001037125
ORF Size:	831 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC214684).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001037125.1</a>
RefSeq Size:	1401 bp
RefSeq ORF:	834 bp
Locus ID:	64718
UniProt ID:	<a href="#">Q9H9P5</a>
Cytogenetics:	16p13.3
Protein Families:	Druggable Genome


[View online »](#)

**MW:** 31.6 kDa

**Gene Summary:** This gene encodes a RING finger protein that may function in Rac signaling. It can bind to Brg/Brm-associated factor 60b and can promote its ubiquitination. Alternative splicing of this gene results in multiple transcript variants. [provided by RefSeq, Jun 2013]