

## Product datasheet for **RC214695A1V**

### Human LINC02897 (NM\_207401) AAV Particle

#### Product data:

|                           |                                                                             |
|---------------------------|-----------------------------------------------------------------------------|
| Product Type:             | AAV Particles                                                               |
| Product Name:             | Human LINC02897 (NM_207401) AAV Particle                                    |
| Tag:                      | Myc-DDK                                                                     |
| Symbol:                   | LINC02897                                                                   |
| Synonyms:                 | FLJ45717                                                                    |
| Mammalian Cell Selection: | None                                                                        |
| Vector:                   | pAAV-AC-Myc-DDK (PS100089)                                                  |
| ORF Nucleotide Sequence:  | >RC214695 representing NM_207401<br>Red=Cloning site Blue=ORF Green=Tags(s) |

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGGCCTGTGCACCCTCCAGCCCTGGGCCCTCCTCGAAATCCTCCACATCCTGCGGGACCTGGACCG  
CCAGCGGCCTCCCAAGCCTGGGCCATTTGCCCCGTCGCCTCCGCCTGGCATTGACTTCTTAGAACCGCG  
GGCAGCAGGGCAGCGTGGGGCAGTGGTGGCTGTCACTCCACGCGGGCCGCCGAGAGGACGCGCCACCC  
CACCCTCCAAACGCCGTCCTCCTGCTTCAGCCAGAGCCGTCCTGACATGGGCCCAGGGACGCGGCTGCA  
GAGGCCACTTCCGGTCACTCCCCGCTGCGGCCTCGCGTCCGGGAGCCGACCCCTCCGCTGCGCCAGCTC  
CGACCGCAGTCTTCGGGAGCAGAAGCAGCGCCGCGCGGGCCCTGACCCACACCCCTCCCCGGCGCCTCCG  
CCTGCAGGACCCCGGCCCTCCCCGGGCTCGCTGGGTCCCTCCGCCCCGGCCGCTCCCGCACCGCGCGCG  
GTGCCTACGAGCTCCAAGGCGGCGCCAGCCAAGACGGTCCCGGCCAAGCGGCGGTGGGTGCGACCCCCAC  
CACGGGCCCCGGCACCGCGGGGAGGGCGCGCTCCTAGGCTGCGGGAGCGGACGGACCCCTCCGACTTCG  
GCCACCTGGCGCCGCGCCTCTCCCCGCCGAAGTCCCTCCCGGGGAGCGCGCCGAACTTTCCGGAGC  
GGGAGCGACTG

**ACGCGT**ACGCGGCGGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

|                          |                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
|--------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Protein Sequence:</b> | <p>&gt;RC214695 representing NM_207401</p> <p>Red=Cloning site Green=Tags(s)</p> <p>MGLCTLQPLGPPRKSSTSCGTWTASGLPSLGHLPRRLRLAFDFLEPRARGQRGGSGGCQSTRAAERTRPP<br/> HPPNAVLLLLQPEPSLTWAQGRGCRGHFRSLPAAASRSGSRTLRCASSDRSLREQKQRRAGPDPTPSPAPP<br/> PAGPRPSPGSLGPSAAPRTARGAYELQGGASQDGPQAAVGATPTTGPGTGGEGALLGCGSGRTPPTS<br/> ATWRRRLLPAEVPPGAAAAANFPERERL</p> <p>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</p>                                                                               |
| <b>Species:</b>          | Human                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>Serotype:</b>         | AAV-2                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>ACCN:</b>             | NM_207401                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>ORF Size:</b>         | 711 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| <b>Buffer:</b>           | PBS with 0.001% Pluronic F68                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| <b>Stability:</b>        | AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus. |
| <b>RefSeq:</b>           | <u>NM_207401.1, NP_997284.1</u>                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>RefSeq Size:</b>      | 2258 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| <b>RefSeq ORF:</b>       | 714 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| <b>Locus ID:</b>         | 388759                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| <b>Cytogenetics:</b>     | 1q44                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| <b>MW:</b>               | 24.2 kDa                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |