

## Product datasheet for **RC214697L2V**

### TRIM14 (NM\_014788) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	TRIM14 (NM_014788) Human Tagged ORF Clone Lentiviral Particle
Symbol:	TRIM14
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_014788
ORF Size:	1326 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC214697).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_014788.2</a>
RefSeq Size:	4454 bp
RefSeq ORF:	1329 bp
Locus ID:	9830
UniProt ID:	<a href="#">Q14142</a>
Cytogenetics:	9q22.33
Domains:	zf-B_box, SPRY, PRY
Protein Families:	Druggable Genome
MW:	49.6 kDa



[View online »](#)

**Gene Summary:**

The protein encoded by this gene is a member of the tripartite motif (TRIM) family. The TRIM motif includes three zinc-binding domains, a RING, a B-box type 1 and a B-box type 2, and a coiled-coil region. The protein localizes to cytoplasmic bodies and its function has not been determined. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Mar 2010]