

Product datasheet for **RC214718A1V**

Human HSPB2 (NM_001541) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human HSPB2 (NM_001541) AAV Particle
Tag: Myc-DDK
Symbol: HSPB2
Synonyms: Hs.78846; HSP27; LOH11CR1K; MKBP
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC214718 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGTCGGGCGCTCAGTGCCACATGCCACCCGGCCACCGCCGAGTACGAATTTGCCAACCCGAGCCGCC
TGGGTGAGCAGCGCTTCGGAGAAGGCCTCCTGCCAGAAGAGATCCTGACCCCACTCTACCATGGCTA
CTATGTCCGGCTCGGGCCGCCAGCTGGGGAGGGCAGCAGGGCAGGGGCTCCGAGCTTAGGCTCAGT
GAGGGCAAGTTCCAGGCATTTCTGGATGTGAGCCACTTACCCAGACGAGGTGACTGTGAGGACTGTGG
ATAACCTGCTGGAGGTGTCTGCCCGGCACCCCGAGCGCCTGGACCGCCACGGCTTCGTGTCCCGAGAGTT
CTGCCGACCTATGTCTGCCTGCTGATGTGACCCCTGGCGAGTCCGAGCTGCTCTCCCATGATGGC
ATCTTAAACCTGGAAGCACCTCGGGGTGGCCGACATTTGGACACAGAGGTCAATGAGGTCTACATCTCCC
TGCTCCCTGCGCTCCTGATCCAGAGGAAGAGGAGGAGGCAGCCATAGTTGAGCCC

ACGCGTACGCGGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC214718 protein sequence
Red=Cloning site Green=Tags(s)

MSGRSVPHAHAPATAEYEFANPSRLGEQRFEGELLPEEILTPTLYHGYVVRPRAAPAGEGSRAGASELRSL
EGKFQAFLDVSHFTPDEVTVRTVDNLLVVSARHPQRLDRHGFVSREFCRTYVLPADVDPWRVRAALSHDG
ILNLEAPRGRHLDTEVNEVYISLLPAPPDPEEEEEAAIVEP

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human
Serotype: AAV-2



[View online »](#)

ACCN:	NM_001541
ORF Size:	546 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001541.2</u>
RefSeq Size:	858 bp
RefSeq ORF:	549 bp
Locus ID:	3316
UniProt ID:	<u>Q16082</u>
Cytogenetics:	11q23.1
MW:	20.2 kDa