

Product datasheet for **RC214772L3V**

CABYR (NM_012189) Human Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | CABYR (NM_012189) Human Tagged ORF Clone Lentiviral Particle |
| Symbol: | CABYR |
| Synonyms: | CABYRa; CABYRc; CABYRc/d; CABYRe; CBP86; CT88; FSP-2; FSP2 |
| Mammalian Cell Selection: | Puromycin |
| Vector: | pLenti-C-Myc-DDK-P2A-Puro (PS100092) |
| Tag: | Myc-DDK |
| ACCN: | NM_012189 |
| ORF Size: | 1479 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RC214772). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_012189.2 |
| RefSeq Size: | 2333 bp |
| RefSeq ORF: | 1482 bp |
| Locus ID: | 26256 |
| UniProt ID: | O75952 |
| Cytogenetics: | 18q11.2 |
| Domains: | R11a |
| MW: | 52.6 kDa |



[View online »](#)

Gene Summary:

To reach fertilization competence, spermatozoa undergo a series of morphological and molecular maturational processes, termed capacitation, involving protein tyrosine phosphorylation and increased intracellular calcium. The protein encoded by this gene localizes to the principal piece of the sperm flagellum in association with the fibrous sheath and exhibits calcium-binding when phosphorylated during capacitation. A pseudogene on chromosome 3 has been identified for this gene. Alternatively spliced transcript variants encoding distinct protein isoforms have been found for this gene. [provided by RefSeq, Jul 2013]