

## Product datasheet for **RC214777L2V**

### **RBM15 (NM\_022768) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	RBM15 (NM_022768) Human Tagged ORF Clone Lentiviral Particle
Symbol:	RBM15
Synonyms:	OTT; OTT1; SPEN
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_022768
ORF Size:	2931 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC214777).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_022768.4</a> , <a href="#">NP_073605.4</a>
RefSeq Size:	3450 bp
RefSeq ORF:	2934 bp
Locus ID:	64783
UniProt ID:	<a href="#">Q96T37</a>
Cytogenetics:	1p13.3
Protein Families:	Transcription Factors
MW:	107 kDa



[View online »](#)

**Gene Summary:**

Members of the SPEN (Split-end) family of proteins, including RBM15, have repressor function in several signaling pathways and may bind to RNA through interaction with spliceosome components (Hiriart et al., 2005 [PubMed 16129689]).[supplied by OMIM, Feb 2009]