

Product datasheet for **RC214826A1V**

Human TXNDC6 (NME9) (NM_178130) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human TXNDC6 (NME9) (NM_178130) AAV Particle
Tag:	Myc-DDK
Symbol:	TXNDC6
Synonyms:	NM23-H9; NXL2; TXL-2; TXL2; TXNDC6
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC214826 representing NM_178130 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGCTCAGTTCCAAAGGACTAACTGTTGTTGATGTCTATCAAGGCTGGTGTGGCCCTGCAAACCTGTGG
TGAGCCTCTTCCAGAAGATGAGGATCGAGGTCGGCCTGGACCTTCTGCACTTTGCATTAGCAGAGGCAGA
TCGTCTTGATGTCCTCGAAAAGTACAGAGGGAAGTGCAGGCCAACCTTTCTGTTTTATGCAATTAAGAT
GAGGCTCTTTCTGATGAAGATGAATGTGTTTCCCATGGAAGAATAATGGTGAAGATGAGGACATGGTTT
CATCAGAGAGGACCTGTACCTTGGCCATCATTAAACCAGATGCAGTGGCCATGGAAGACTGATGAGAT
TATCATGAAGATTACAGGAAGCTGGGTTTGAAATTCTAACAATGAAGAGAGAACCATGACAGAGGCAGAA
GTGCGACTTTTCTACCAACACAAAGCTGGAGAGGAGGCATTTGAGAAGCTGGTACATCACATGTGCAGTG
GACCAAGCCACCTCCTGATCCTCACCAGGACTGAGGGCTTCGAGGACGTGGTCACTACCTGGCGAACCGT
CATGGGCCCCGTGACCCCAATGTGGCCAGGAGGGAGCAGCCAGAAAGTCTCCGAGCTCAGTACGGCACA
GAAATGCCCTTCAATGCCGTCCATGGAAGCCGGGACAGAGAAGATGCTGACAGAGAACTGGCATTGCTCT
TCCCCAGTTTGAAATTTTCAGACAAAGATACAGAAGCCCTCAGGGTGAGTCCAGTACACAGCCAAGTT
GAAGATCACTGACCTAGAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC214826 representing NM_178130 Red=Cloning site Green=Tags(s) MLSSKGLTVVDVYQGWCGPCKPVVSLFQKMRIEVLDLLHFALEADRLDVLEKYRGKCEPTFLFYAIKD EALSDEDECVSHGKNNGEDEDMVSSERTCTLAIIKPDVAHGKTDEIIMKIQEAGFEILTNEERTMTEAE VRLFYQHKAGEEAFEKLVHHMCSGSPSHLLILTRTEGFEDVTTWRTVMGPRDPNVARREQPESLRAQYGT EMPFNAVHGSRDREDADRELALLFPSLKFSDKDTEAPQGESSTQPRLKITDLD TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_178130
ORF Size:	789 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_178130.2</u>
RefSeq Size:	2133 bp
RefSeq ORF:	792 bp
Locus ID:	347736
UniProt ID:	<u>Q86XW9</u>
Cytogenetics:	3q22.3
MW:	29.5 kDa