

## Product datasheet for **RC214828L4V**

### ZNF266 (NM\_006631) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	ZNF266 (NM_006631) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ZNF266
Synonyms:	HZF1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_006631
ORF Size:	1647 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC214828).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_006631.2</a>
RefSeq Size:	3790 bp
RefSeq ORF:	1650 bp
Locus ID:	10781
UniProt ID:	<a href="#">Q14584</a>
Cytogenetics:	19p13.2
Domains:	zf-C2H2
Protein Families:	Transcription Factors



[View online »](#)

**MW:** 62.1 kDa

**Gene Summary:** This gene encodes a protein containing many tandem zinc-finger motifs. Zinc fingers are protein or nucleic acid-binding domains, and may be involved in a variety of functions, including regulation of transcription. This gene is located in a cluster of similar genes encoding zinc finger proteins on chromosome 19. Alternative splicing results in multiple transcript variants for this gene. [provided by RefSeq, Sep 2012]