

Product datasheet for **RC214846L4V**

TMEM132A (NM_017870) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	TMEM132A (NM_017870) Human Tagged ORF Clone Lentiviral Particle
Symbol:	TMEM132A
Synonyms:	GBP; HSPA5BP1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_017870
ORF Size:	3072 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC214846).
OTI Disclaimer:	<p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info</p>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_017870.2
RefSeq Size:	3486 bp
RefSeq ORF:	3075 bp



[View online »](#)

Locus ID:	54972
UniProt ID:	Q24JP5
Cytogenetics:	11q12.2
Protein Families:	Transmembrane
MW:	110 kDa
Gene Summary:	This gene encodes a protein that is highly similar to the rat Grp78-binding protein (GBP). Alternatively spliced transcript variants encoding different isoforms have been described. [provided by RefSeq, Jul 2008]