

## Product datasheet for **RC214861L2V**

### SPINT1 (NM\_003710) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	SPINT1 (NM_003710) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SPINT1
Synonyms:	HAI; HAI1; MANSC2
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_003710
ORF Size:	1539 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC214861).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_003710.3</a>
RefSeq Size:	2436 bp
RefSeq ORF:	1542 bp
Locus ID:	6692
UniProt ID:	<a href="#">O43278</a>
Cytogenetics:	15q15.1
Domains:	Idl_recept_a, KU
Protein Families:	Druggable Genome, Secreted Protein, Transmembrane



[View online »](#)

**MW:** 56.89 kDa

**Gene Summary:** The protein encoded by this gene is a member of the Kunitz family of serine protease inhibitors. The protein is a potent inhibitor specific for HGF activator and is thought to be involved in the regulation of the proteolytic activation of HGF in injured tissues. Alternative splicing results in multiple variants encoding different isoforms. [provided by RefSeq, Jul 2008]