

Product datasheet for **RC214873L2V**

CMA1 (NM_001836) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	CMA1 (NM_001836) Human Tagged ORF Clone Lentiviral Particle
Symbol:	CMA1
Synonyms:	chymase; CYH; MCT1
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_001836
ORF Size:	741 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC214873).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001836.1
RefSeq Size:	744 bp
RefSeq ORF:	744 bp
Locus ID:	1215
UniProt ID:	P23946
Cytogenetics:	14q12
Protein Families:	Druggable Genome, Protease, Secreted Protein
Protein Pathways:	Renin-angiotensin system



[View online »](#)

MW: 27.1 kDa

Gene Summary: This gene encodes a chymotryptic serine proteinase that belongs to the peptidase family S1. It is expressed in mast cells and is thought to function in the degradation of the extracellular matrix, the regulation of submucosal gland secretion, and the generation of vasoactive peptides. In the heart and blood vessels, this protein, rather than angiotensin converting enzyme, is largely responsible for converting angiotensin I to the vasoactive peptide angiotensin II. Alternative splicing results in multiple variants. [provided by RefSeq, Apr 2015]