

Product datasheet for **RC214903L4V**

Adenosine A3 Receptor (ADORA3) (NM_001081976) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Adenosine A3 Receptor (ADORA3) (NM_001081976) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Adenosine A3 Receptor
Synonyms:	AD026
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001081976
ORF Size:	798 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC214903).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001081976.1 , NP_001075445.1
RefSeq Size:	1185 bp
RefSeq ORF:	801 bp
Locus ID:	57413
UniProt ID:	P33765
Cytogenetics:	1p13.2
MW:	30.1 kDa



[View online »](#)

Gene Summary:

This gene encodes a transmembrane and immunoglobulin domain-containing protein. Alternative splicing results in multiple transcript variants, one of which shares its 5' terminal exon with that of the overlapping adenosine A3 receptor gene (GenelD:140), thus resulting in a fusion product. [provided by RefSeq, Nov 2014]