

## Product datasheet for **RC214916L4V**

### ATP5SL (DMAC2) (NM\_018035) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	ATP5SL (DMAC2) (NM_018035) Human Tagged ORF Clone Lentiviral Particle
Symbol:	DMAC2
Synonyms:	ATP5SL
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_018035
ORF Size:	771 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC214916).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_018035.1</a> , <a href="#">NP_060505.1</a>
RefSeq Size:	1692 bp
RefSeq ORF:	774 bp
Locus ID:	55101
UniProt ID:	<a href="#">Q9NW81</a>
Cytogenetics:	19q13.2
MW:	29.1 kDa



[View online »](#)

**Gene Summary:**

Required for the assembly of the mitochondrial NADH:ubiquinone oxidoreductase complex (complex I). Involved in the assembly of the distal region of complex I.[UniProtKB/Swiss-Prot Function]