

## Product datasheet for **RC214943L3V**

### TAC1 (NM\_013996) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	TAC1 (NM_013996) Human Tagged ORF Clone Lentiviral Particle
Symbol:	TAC1
Synonyms:	Hs.2563; NK2; NKNA; NPK; TAC2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_013996
ORF Size:	333 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC214943).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_013996.1</a>
RefSeq Size:	1048 bp
RefSeq ORF:	336 bp
Locus ID:	6863
UniProt ID:	<a href="#">P20366</a>
Cytogenetics:	7q21.3
Domains:	Tachykinin
Protein Families:	Druggable Genome, Secreted Protein



[View online »](#)

**MW:** 13.04 kDa

**Gene Summary:** This gene encodes four products of the tachykinin peptide hormone family, substance P and neurokinin A, as well as the related peptides, neuropeptide K and neuropeptide gamma. These hormones are thought to function as neurotransmitters which interact with nerve receptors and smooth muscle cells. They are known to induce behavioral responses and function as vasodilators and secretagogues. Substance P is an antimicrobial peptide with antibacterial and antifungal properties. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Nov 2014]