

Product datasheet for **RC214965L4V**

ZNF264 (NM_003417) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	ZNF264 (NM_003417) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ZNF264
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_003417
ORF Size:	1881 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC214965).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_003417.2
RefSeq Size:	12172 bp
RefSeq ORF:	1884 bp
Locus ID:	9422
UniProt ID:	O43296
Cytogenetics:	19q13.43
Domains:	KRAB, zf-C2H2
Protein Families:	Transcription Factors
MW:	70.5 kDa



[View online »](#)

Gene Summary:

This gene encodes a zinc finger protein and belongs to the krueppel C2H2-type zinc-finger protein family. Zinc finger proteins are often localized in the nucleus, bind nucleic acids, and regulate transcription. [provided by RefSeq, Jan 2010]