

Product datasheet for **RC215000L3V**

BATF3 (NM_018664) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	BATF3 (NM_018664) Human Tagged ORF Clone Lentiviral Particle
Symbol:	BATF3
Synonyms:	JDP1; JUNDM1; SNFT
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_018664
ORF Size:	381 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC215000).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_018664.1
RefSeq Size:	992 bp
RefSeq ORF:	384 bp
Locus ID:	55509
UniProt ID:	Q9NR55
Cytogenetics:	1q32.3
MW:	14.5 kDa



[View online »](#)

Gene Summary:

This gene encodes a member of the basic leucine zipper protein family. The encoded protein functions as a transcriptional repressor when heterodimerizing with JUN. The protein may play a role in repression of interleukin-2 and matrix metalloproteinase-1 transcription. [provided by RefSeq, Feb 2009]