

## Product datasheet for **RC215022L4V**

### ERMAP (NM\_001017922) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	ERMAP (NM_001017922) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ERMAP
Synonyms:	BTN5; PRO2801; RD; SC
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001017922
ORF Size:	1425 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC215022).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001017922.1</a>
RefSeq Size:	3485 bp
RefSeq ORF:	1428 bp
Locus ID:	114625
UniProt ID:	<a href="#">Q96PL5</a>
Cytogenetics:	1p34.2
Protein Families:	Druggable Genome, Transmembrane
MW:	52.6 kDa



[View online »](#)

**Gene Summary:**

The protein encoded by this gene is a cell surface transmembrane protein that may act as an erythroid cell receptor, possibly as a mediator of cell adhesion. Polymorphisms in this gene are responsible for the Scianna/Radin blood group system. Two transcript variants encoding the same protein have been found for this gene. [provided by RefSeq, Jul 2008]