

## Product datasheet for **RC215028L4V**

### LYK5 (STRADA) (NM\_001003788) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

|                           |  |
|---------------------------|--|
| Product Type:             | Lentiviral Particles   |
| Product Name:             | LYK5 (STRADA) (NM_001003788) Human Tagged ORF Clone Lentiviral Particle  |
| Symbol:                   | LYK5   |
| Synonyms:                 | LYK5; NY-BR-96; PMSE; StIk; STRAD; STRAD alpha   |
| Mammalian Cell Selection: | Puromycin  |
| Vector:                   | pLenti-C-mGFP-P2A-Puro (PS100093)  |
| Tag:                      | mGFP   |
| ACCN:                     | NM_001003788   |
| ORF Size:                 | 1119 bp  |
| ORF Nucleotide Sequence:  | The ORF insert of this clone is exactly the same as(RC215028).   |
| OTI Disclaimer:           | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a> |
| OTI Annotation:           | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.   |
| RefSeq:                   | <a href="#">NM_001003788.1</a>   |
| RefSeq Size:              | 2194 bp  |
| RefSeq ORF:               | 1122 bp  |
| Locus ID:                 | 92335  |
| UniProt ID:               | <a href="#">Q7RTN6</a>   |
| Cytogenetics:             | 17q23.3  |
| Protein Families:         | Druggable Genome, Protein Kinase   |
| Protein Pathways:         | mTOR signaling pathway   |



[View online »](#)

**MW:** 41.5 kDa

**Gene Summary:** The protein encoded by this gene contains a STE20-like kinase domain, but lacks several residues that are critical for catalytic activity, so it is termed a 'pseudokinase'. The protein forms a heterotrimeric complex with serine/threonine kinase 11 (STK11, also known as LKB1) and the scaffolding protein calcium binding protein 39 (CAB39, also known as MO25). The protein activates STK11 leading to the phosphorylation of both proteins and excluding STK11 from the nucleus. The protein is necessary for STK11-induced G1 cell cycle arrest. A mutation in this gene has been shown to result in polyhydramnios, megalencephaly, and symptomatic epilepsy (PMSE) syndrome. Multiple transcript variants encoding different isoforms have been found for this gene. Additional transcript variants have been described but their full-length nature is not known. [provided by RefSeq, Sep 2009]