

Product datasheet for RC215069L3

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Complement C3 (C3) (NM_000064) Human Tagged Lenti ORF Clone

Product data:

Product Type: Expression Plasmids

Product Name: Complement C3 (C3) (NM_000064) Human Tagged Lenti ORF Clone

Tag: Myc-DDK

Symbol: Complement C3

Synonyms: AHUS5; ARMD9; ASP; C3a; C3b; CPAMD1; HEL-S-62p

Mammalian Cell Puromycin

Selection:

Vector: pLenti-C-Myc-DDK-P2A-Puro (PS100092)

E. coli Selection: Chloramphenicol (34 ug/mL)

ORF Nucleotide The ORF insert of this clone

Sequence:

The ORF insert of this clone is exactly the same as(RC215069).

Restriction Sites: Sgfl-Mlul

Cloning Scheme:





^{*} The last codon before the Stop codon of the ORF.

ACCN: NM_000064

ORF Size: 4989 bp





Complement C3 (C3) (NM_000064) Human Tagged Lenti ORF Clone - RC215069L3

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of

reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube

containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method: 1. Centrifuge at 5,000xg for 5min.

2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.

3. Close the tube and incubate for 10 minutes at room temperature.

4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid

at the bottom.

5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of

shipping when stored at -20°C.

RefSeq: <u>NM 000064.1</u>, <u>NP 000055.1</u>

RefSeq Size: 5067 bp
RefSeq ORF: 4992 bp
Locus ID: 718

 UniProt ID:
 P01024

 Cytogenetics:
 19p13.3

Domains: ANATO, NTR, A2M, A2M N

Protein Families: Druggable Genome

Protein Pathways: Complement and coagulation cascades, Systemic lupus erythematosus

MW: 187 kDa


Gene Summary: Complement component C3 plays a central role in the activation of complement system. Its

activation is required for both classical and alternative complement activation pathways. The encoded preproprotein is proteolytically processed to generate alpha and beta subunits that form the mature protein, which is then further processed to generate numerous peptide products. The C3a peptide, also known as the C3a anaphylatoxin, modulates inflammation and possesses antimicrobial activity. Mutations in this gene are associated with atypical hemolytic uremic syndrome and age-related macular degeneration in human patients.


[provided by RefSeq, Nov 2015]



Product images:



Circular map for RC215069L3



Double digestion of RC215069L3 using Sgfl and Mlul