

## **Product datasheet for RC215074L1**

## OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

### Amylin (IAPP) (NM\_000415) Human Tagged Lenti ORF Clone

#### **Product data:**

**Product Type:** Expression Plasmids

**Product Name:** Amylin (IAPP) (NM\_000415) Human Tagged Lenti ORF Clone

Tag: Myc-DDK
Symbol: Amylin
Synonyms: DAP; IAP

Mammalian Cell None

Selection:

Vector:pLenti-C-Myc-DDK (PS100064)E. coli Selection:Chloramphenicol (34 ug/mL)

ORF Nucleotide The ORF insert of this clone is exactly the same as(RC215074).

Sequence:

**Restriction Sites:** Sgfl-Mlul

**Cloning Scheme:** 





<sup>\*</sup> The last codon before the Stop codon of the ORF.

**ACCN:** NM\_000415

ORF Size: 267 bp





### Amylin (IAPP) (NM\_000415) Human Tagged Lenti ORF Clone - RC215074L1

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of

reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube

containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:** 1. Centrifuge at 5,000xg for 5min.

2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.

3. Close the tube and incubate for 10 minutes at room temperature.

4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid

at the bottom.

5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of

shipping when stored at -20°C.

**RefSeq:** <u>NM 000415.1</u>

 RefSeq Size:
 1462 bp

 RefSeq ORF:
 270 bp

 Locus ID:
 3375

 UniProt ID:
 P10997

Cytogenetics: 12p12.1

Protein Families: Druggable Genome, Secreted Protein
Protein Pathways: Maturity onset diabetes of the young

**MW:** 9.81 kDa


**Gene Summary:** This gene encodes a member of the calcitonin family of peptide hormones. This hormone is

released from pancreatic beta cells following food intake to regulate blood glucose levels and act as a satiation signal. Human patients with type 1 and advanced type 2 diabetes exhibit reduced levels of the encoded hormone in blood and pancreas. This protein also exhibits a


bactericidal, antimicrobial activity. [provided by RefSeq, Jul 2016]



# **Product images:**



Circular map for RC215074L1



Double digestion of RC215074L1 using Sgfl and Mlul