

Product datasheet for **RC215080L4V**

ZNF366 (NM_152625) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	ZNF366 (NM_152625) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ZNF366
Synonyms:	DC-SCRIPT; DCSCRIPT
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_152625
ORF Size:	2232 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC215080).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_152625.1
RefSeq Size:	2770 bp
RefSeq ORF:	2235 bp
Locus ID:	167465
UniProt ID:	Q8N895
Cytogenetics:	5q13.2
Domains:	zf-C2H2
Protein Families:	Transcription Factors


[View online »](#)

MW: 85.1 kDa

Gene Summary: Has transcriptional repression activity. Acts as corepressor of ESR1; the function seems to involve CTBP1 and histone deacetylases.[UniProtKB/Swiss-Prot Function]