

Product datasheet for **RC215092L3V**

ICAM4 (NM_022377) Human Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | ICAM4 (NM_022377) Human Tagged ORF Clone Lentiviral Particle |
| Symbol: | ICAM4 |
| Synonyms: | CD242; LW |
| Mammalian Cell Selection: | Puromycin |
| Vector: | pLenti-C-Myc-DDK-P2A-Puro (PS100092) |
| Tag: | Myc-DDK |
| ACCN: | NM_022377 |
| ORF Size: | 711 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RC215092). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_022377.1 , NP_071772.1 |
| RefSeq Size: | 1390 bp |
| RefSeq ORF: | 714 bp |
| Locus ID: | 3386 |
| Cytogenetics: | 19p13.2 |
| Protein Families: | Secreted Protein, Transmembrane |
| MW: | 23.6 kDa |



[View online »](#)

Gene Summary:

This gene encodes the Landsteiner-Wiener (LW) blood group antigen(s) that belongs to the immunoglobulin (Ig) superfamily, and that shares similarity with the intercellular adhesion molecule (ICAM) protein family. This ICAM protein contains 2 Ig-like C2-type domains and binds to the leukocyte adhesion LFA-1 protein. The molecular basis of the LW(A)/LW(B) blood group antigens is a single aa variation at position 100; Gln-100=LW(A) and Arg-100=LW(B). Alternative splicing results in multiple transcript variants encoding distinct isoforms. [provided by RefSeq, Jul 2008]