

Product datasheet for **RC215098L2V**

Striatin (STRN) (NM_003162) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Striatin (STRN) (NM_003162) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Striatin
Synonyms:	PPP2R6A; SG2NA; STRN1
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_003162
ORF Size:	2340 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC215098).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_003162.2
RefSeq Size:	3479 bp
RefSeq ORF:	2343 bp
Locus ID:	6801
UniProt ID:	O43815
Cytogenetics:	2p22.2
Domains:	WD40
MW:	86 kDa



[View online »](#)

Gene Summary:

Calmodulin-binding protein which may function as scaffolding or signaling protein and may play a role in dendritic Ca(2+) signaling.[UniProtKB/Swiss-Prot Function]