

Product datasheet for **RC215146L4V**

APBA1 (NM_001163) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	APBA1 (NM_001163) Human Tagged ORF Clone Lentiviral Particle
Symbol:	APBA1
Synonyms:	D9S411E; LIN10; MINT1; X11; X11A; X11ALPHA
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001163
ORF Size:	2511 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC215146).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001163.3
RefSeq Size:	6603 bp
RefSeq ORF:	2514 bp
Locus ID:	320
UniProt ID:	Q02410
Cytogenetics:	9q21.12
Domains:	PDZ, PID
MW:	92.9 kDa



[View online »](#)

Gene Summary:

The protein encoded by this gene is a member of the X11 protein family. It is a neuronal adapter protein that interacts with the Alzheimer's disease amyloid precursor protein (APP). It stabilizes APP and inhibits production of proteolytic APP fragments including the A beta peptide that is deposited in the brains of Alzheimer's disease patients. This gene product is believed to be involved in signal transduction processes. It is also regarded as a putative vesicular trafficking protein in the brain that can form a complex with the potential to couple synaptic vesicle exocytosis to neuronal cell adhesion. [provided by RefSeq, Jul 2008]