

## Product datasheet for **RC215162L3V**

### TRAF7 (NM\_032271) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	TRAF7 (NM_032271) Human Tagged ORF Clone Lentiviral Particle
Symbol:	TRAF7
Synonyms:	CAFADADD; RFWD1; RNF119
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_032271
ORF Size:	2010 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC215162).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_032271.2</a>
RefSeq Size:	3676 bp
RefSeq ORF:	2013 bp
Locus ID:	84231
UniProt ID:	<a href="#">Q6Q0C0</a>
Cytogenetics:	16p13.3
Domains:	zf-TRAF, WD40, RING
Protein Families:	Druggable Genome, Transcription Factors



[View online »](#)

**MW:** 74.4 kDa

**Gene Summary:** Tumor necrosis factor (TNF; see MIM 191160) receptor-associated factors, such as TRAF7, are signal transducers for members of the TNF receptor superfamily (see MIM 191190). TRAFs are composed of an N-terminal cysteine/histidine-rich region containing zinc RING and/or zinc finger motifs; a coiled-coil (leucine zipper) motif; and a homologous region that defines the TRAF family, the TRAF domain, which is involved in self-association and receptor binding. [supplied by OMIM, Apr 2004]