

## Product datasheet for **RC215172L2V**

### HNT (NTM) (NM\_016522) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	HNT (NTM) (NM_016522) Human Tagged ORF Clone Lentiviral Particle
Symbol:	HNT
Synonyms:	CEPU-1; HNT; IGLON2; NTRI
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_016522
ORF Size:	1032 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC215172).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_016522.2</a> , <a href="#">NP_057606.1</a>
RefSeq Size:	3381 bp
RefSeq ORF:	1035 bp
Locus ID:	50863
UniProt ID:	<a href="#">Q9P121</a>
Cytogenetics:	11q25
Domains:	ig
Protein Families:	Transmembrane



[View online »](#)

**MW:** 37.8 kDa

**Gene Summary:** This gene encodes a member of the IgLON (LAMP, OBCAM, Ntm) family of immunoglobulin (Ig) domain-containing glycosylphosphatidylinositol (GPI)-anchored cell adhesion molecules. The encoded protein may promote neurite outgrowth and adhesion via a homophilic mechanism. This gene is closely linked to a related family member, opioid binding protein/cell adhesion molecule-like (OPCML), on chromosome 11. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jan 2009]