

Product datasheet for **RC215177L4V**

CABLES2 (NM_031215) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	CABLES2 (NM_031215) Human Tagged ORF Clone Lentiviral Particle
Symbol:	CABLES2
Synonyms:	C20orf150; dj908M14.2; ik3-2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_031215
ORF Size:	1434 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC215177).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_031215.2 , NP_112492.2
RefSeq Size:	3785 bp
RefSeq ORF:	1437 bp
Locus ID:	81928
UniProt ID:	Q9BTV7
Cytogenetics:	20q13.33
MW:	52.1 kDa
Gene Summary:	Unknown. Probably involved in G1-S cell cycle transition.[UniProtKB/Swiss-Prot Function]


[View online »](#)