

## Product datasheet for **RC215215L4V**

### **PCK2 (NM\_001018073) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	PCK2 (NM_001018073) Human Tagged ORF Clone Lentiviral Particle
Symbol:	PCK2
Synonyms:	PEPCK; PEPCK-M; PEPCK2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001018073
ORF Size:	1323 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC215215).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001018073.1</a> , <a href="#">NP_001018083.1</a>
RefSeq Size:	1730 bp
RefSeq ORF:	1326 bp
Locus ID:	5106
UniProt ID:	<a href="#">Q16822</a>
Cytogenetics:	14q11.2-q12
Protein Families:	ES Cell Differentiation/IPS



[View online »](#)

<b>Protein Pathways:</b>	Adipocytokine signaling pathway, Citrate cycle (TCA cycle), Glycolysis / Gluconeogenesis, Insulin signaling pathway, Metabolic pathways, PPAR signaling pathway, Pyruvate metabolism
<b>MW:</b>	47.56 kDa
<b>Gene Summary:</b>	This gene encodes a mitochondrial enzyme that catalyzes the conversion of oxaloacetate to phosphoenolpyruvate in the presence of guanosine triphosphate (GTP). A cytosolic form of this protein is encoded by a different gene and is the key enzyme of gluconeogenesis in the liver. Alternatively spliced transcript variants have been described. [provided by RefSeq, Apr 2014]