

Product datasheet for **RC215222L1V**

EBF3 (NM_001005463) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	EBF3 (NM_001005463) Human Tagged ORF Clone Lentiviral Particle
Symbol:	EBF3
Synonyms:	COE3; EBF-3; HADDS; O/E-2; OE-2
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_001005463
ORF Size:	1653 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC215222).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001005463.1
RefSeq Size:	4412 bp
RefSeq ORF:	1656 bp
Locus ID:	253738
UniProt ID:	Q9H4W6
Cytogenetics:	10q26.3
Protein Families:	Transcription Factors
MW:	60.4 kDa



[View online »](#)

Gene Summary:

This gene encodes a member of the early B-cell factor (EBF) family of DNA binding transcription factors. EBF proteins are involved in B-cell differentiation, bone development and neurogenesis, and may also function as tumor suppressors. The encoded protein inhibits cell survival through the regulation of genes involved in cell cycle arrest and apoptosis, and aberrant methylation or deletion of this gene may play a role in multiple malignancies including glioblastoma multiforme and gastric carcinoma. [provided by RefSeq, Sep 2011]