

## Product datasheet for **RC215223L2V**

### FECH (NM\_001012515) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	FECH (NM_001012515) Human Tagged ORF Clone Lentiviral Particle
Symbol:	FECH
Synonyms:	EPP; EPP1; FCE
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_001012515
ORF Size:	1287 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC215223).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001012515.1</a>
RefSeq Size:	3835 bp
RefSeq ORF:	1290 bp
Locus ID:	2235
UniProt ID:	<a href="#">P22830</a>
Cytogenetics:	18q21.31
Protein Families:	Druggable Genome
Protein Pathways:	Metabolic pathways, Porphyrin and chlorophyll metabolism



[View online »](#)

**MW:** 48.63 kDa

**Gene Summary:** The protein encoded by this gene is localized to the mitochondrion, where it catalyzes the insertion of the ferrous form of iron into protoporphyrin IX in the heme synthesis pathway. Mutations in this gene are associated with erythropoietic protoporphyria. Two transcript variants encoding different isoforms have been found for this gene. A pseudogene of this gene is found on chromosome 3.[provided by RefSeq, May 2010]