

Product datasheet for **RC215229L4V**

SLC39A4 (NM_017767) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	SLC39A4 (NM_017767) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SLC39A4
Synonyms:	AEZ; AWMS2; ZIP4
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_017767
ORF Size:	1866 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC215229).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_017767.2 , NP_060237.2
RefSeq Size:	2317 bp
RefSeq ORF:	1869 bp
Locus ID:	55630
UniProt ID:	Q6P5W5
Cytogenetics:	8q24.3
Protein Families:	Transmembrane
MW:	66 kDa



[View online »](#)

Gene Summary:

This gene encodes a member of the zinc/iron-regulated transporter-like protein (ZIP) family. The encoded protein localizes to cell membranes and is required for zinc uptake in the intestine. Mutations in this gene result in acrodermatitis enteropathica. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2013]