

Product datasheet for **RC215233L4V**

HDAC7 (NM_015401) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	HDAC7 (NM_015401) Human Tagged ORF Clone Lentiviral Particle
Symbol:	HDAC7
Synonyms:	HD7; HD7A; HDAC7A
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_015401
ORF Size:	2973 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC215233).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_015401.3
RefSeq Size:	4301 bp
RefSeq ORF:	2976 bp
Locus ID:	51564
UniProt ID:	Q8WUI4
Cytogenetics:	12q13.11
Domains:	Hist_deacetyl
Protein Families:	Druggable Genome, Transcription Factors



[View online »](#)

MW: 106.6 kDa

Gene Summary: Histones play a critical role in transcriptional regulation, cell cycle progression, and developmental events. Histone acetylation/deacetylation alters chromosome structure and affects transcription factor access to DNA. The protein encoded by this gene has sequence homology to members of the histone deacetylase family. This gene is orthologous to mouse HDAC7 gene whose protein promotes repression mediated via the transcriptional corepressor SMRT. Alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008]