

Product datasheet for RC215235L4

IL21 (NM_021803) Human Tagged Lenti ORF Clone

Product data:

Product Type: Expression Plasmids

Product Name: IL21 (NM_021803) Human Tagged Lenti ORF Clone

Tag: mGFP Symbol: IL21

Synonyms: CVID11; IL-21; Za11

Mammalian Cell Puromycin

Selection:

Vector: pLenti-C-mGFP-P2A-Puro (PS100093)

E. coli Selection: Chloramphenicol (34 ug/mL)

ORF Nucleotide The ORF insert of this clone is exactly the same as(RC215235).

Sequence:

Restriction Sites: Sgfl-Mlul

Cloning Scheme:





^{*} The last codon before the Stop codon of the ORF

ACCN: NM_021803

ORF Size: 486 bp



OriGene Technologies, Inc. 9620 Medical Center Drive, Ste 200

CN: techsupport@origene.cn

Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com

IL21 (NM_021803) Human Tagged Lenti ORF Clone - RC215235L4

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of

reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube

containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method: 1. Centrifuge at 5,000xg for 5min.

2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.

3. Close the tube and incubate for 10 minutes at room temperature.

4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid

at the bottom.

5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of

shipping when stored at -20°C.

RefSeq: <u>NM 021803.1</u>

RefSeq Size: 642 bp
RefSeq ORF: 489 bp
Locus ID: 59067
UniProt ID: Q9HBE4

Cytogenetics:

Protein Families: Druggable Genome, Secreted Protein

4q27

Protein Pathways: Cytokine-cytokine receptor interaction, Jak-STAT signaling pathway

MW: 18.5 kDa


Gene Summary: This gene encodes a member of the common-gamma chain family of cytokines with

immunoregulatory activity. The encoded protein plays a role in both the innate and adaptive immune responses by inducing the differentiation, proliferation and activity of multiple target cells including macrophages, natural killer cells, B cells and cytotoxic T cells. Dysregulation of this gene plays a role in multiple immune-mediated diseases including lupus, psoriasis and chronic inflammatory diseases. Alternatively spliced transcript variants encoding multiple


isoforms have been observed for this gene. [provided by RefSeq, Dec 2011]



Product images:



Circular map for RC215235L4



Double digestion of RC215235L4 using Sgfl and Mlul