

Product datasheet for **RC215280**

OSCP1 (NM_206837) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	OSCP1 (NM_206837) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	OSCP1
Synonyms:	C1orf102; NOR1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>RC215280 representing NM_206837 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGTCGGTGCGGACGCTACCGCTGCTCTTCTTGAACCTGGGCGGGGAGATGCTTTACATCCTCGACCAAC
GGCTGCGGGCCCAGAACATCCCGGGAGACAAGGCCCGCAAAGTTCTGAATGACATCATCTCCACCATGTT
CAATAGAAAGTTTATGGAGGAATTATCAAGCCTCAAGAGCTCTACTCCAAGAAGCCCTGAGGACTGTC
TATGAGCGCCTGGCTCATGCCTCCATTATGAACTGAACCAGGCCAGCATGGATAAGCTCTATGACCTGA
TGACCATGGCTTTCAAATATCAAGTATTGCTGTGTCCCCGACCAAGGATGTGCTGCTGGTCACTTTCAA
TCACTTGATACCATCAAGGGATTCACTCCGAGACTCCCCAACCATCCTGCAGCAAGTGGACGAGACTTTG
CGGCAGCTGACAGAAATATATGGTGGTCTCTCTGCAGGGGAGTCCAGCTGATCCGGCAGACACTCCTCA
TCTTCTTCCAAGACCTGCACATCCGAAGAGAGTGCTTCTGGAAGCATATGCAAGAACTGTTACAGTGC
ACAGTGTTCCTTCCAGGAGTGGGAACTATATGGAGGACAGGAGAAGGAGACTGACTTTCCATTATGTA
TACCCTTTGAGAACTTTGAGTTTTGTACCTGGGCAAAA

ACGCGTACGCGGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_206837.3](#)

RefSeq Size: 816 bp

RefSeq ORF: 672 bp

Locus ID: 127700


UniProt ID: [Q8WVF1](#)

Cytogenetics: 1p34.3

MW: 26.4 kDa

Gene Summary: May be involved in drug clearance in the placenta.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RC215280