

Product datasheet for **RC215372L3V**

FADS6 (NM_178128) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	FADS6 (NM_178128) Human Tagged ORF Clone Lentiviral Particle
Symbol:	FADS6
Synonyms:	FP18279
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_178128
ORF Size:	1050 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC215372).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_178128.3
RefSeq Size:	2052 bp
RefSeq ORF:	1107 bp
Locus ID:	283985
UniProt ID:	Q8N9I5
Cytogenetics:	17q25.1
Protein Families:	Transmembrane
MW:	39.8 kDa



[View online »](#)